



Equipos automáticos para ERS, ELISA, lectores para contadores hematológicos.

## Heal Force ® Highspeed Refrigerated Benchtop Centrifuge Neofuge 23R

Neofuge 23R was developed under the brand name of Heal Force, able to get excellent centrifugal effect and easy to use with high reliability and safety, it can be widely used in molecular biology and biochemical research, medical diagnosis etc.

Microprocessor control, all parameters are set by control panel

Max. speed 23300rpm, Max force 50377xg

Multi-purpose, 8 rotors to select

rotor can keep in low temperature at Max. speed

9 acceleration/deceleration curves

low-noise brushless motor, maintenance-free



### Specifications

<b>Min/Max. speed</b>	<b>rpm</b>	<b>300~23300, adjustable in steplength of 10</b>
<b>Max. RCF</b>	<b>xg</b>	<b>50377</b>
<b>Max.capacity</b>	<b>ml</b>	<b>4x180</b>
<b>drive</b>		<b>induction drive without carbon brushes</b>
<b>control</b>		<b>Microprocessor control</b>
<b>Braking/acceleration curves</b>		<b>9/9</b>
<b>programs</b>		<b>9</b>
<b>Noise at max. speed</b>	<b>dBA</b>	<b>≤66</b>
<b>Timer setting</b>		<b>≤9h59min</b>
<b>Temperature range</b>		<b>20 40</b>
<b>Rotor recognition</b>		<b>automatic</b>
<b>Tiimbalance detection</b>		<b>automatic</b>
<b>Equipment safety</b>		<b>Lid lock and interlock</b>
<b>Dimesions(h×w×d)</b>	<b>mm</b>	<b>400×765×590</b>
<b>Dimesions crated (h×w×d)</b>	<b>mm</b>	<b>670×980×810</b>
<b>Weight net/gross (excluding rotor)</b>	<b>kg</b>	<b>150/180</b>
<b>Power consumption</b>	<b>w</b>	<b>1800</b>







